


**HeartSine**<sup>®</sup>  
samaritan<sup>®</sup> PAD  
Family

Accessories

**stryker**



**Adult Pad-Pak**  
For patients > 8 years or  
25 kg (55 lb).

**Pad-Pak-03**  
**11516-000003**



**Pediatric-Pak™**  
For patients 1 to 8 years or  
up to 25 kg (55 lb).

**Pad-Pak-04**  
**11516-000004**



**Adult Pad-Pak**  
With TSO-C142a certification  
for use on aircraft.

**Pad-Pak-07**  
**11516-000027**



**USB Data Cable**  
For use with the Saver EVO™  
data management software  
to download event data and to  
perform software upgrades to  
the HeartSine samaritan PAD  
non-connected devices.

**PAD-ACC-02**  
**11516-000018**

## Cabinets and Mounting Options



**Rotaid Plus Wall Cabinet  
with Alarm**  
For indoor use.

**11996-000442** (Green)  
**11996-000441** (White)  
**11996-000443** (Red)  
**11996-000444** (Yellow)




**Rotaid Solid Plus Wall  
Cabinet with Alarm**  
For indoor or outdoor use.

**11996-000446** (Green)  
**11996-000445** (White)  
**11996-000447** (Red)  
**11996-000448** (Yellow)



**Rotaid Solid Plus Heat  
Wall Cabinet with Heat  
and Alarm**  
For outdoor use.

**11996-000450** (Green)  
**11996-000449** (White)  
**11996-000451** (Red)  
**11996-000452** (Yellow)



**AED Wall Bracket**  
(Metal)  
Dimensions:  
15.6 H x 17.5 W x 11.4 D cm  
(6.14 H x 6.88 W x 4.48 D in)

**PAD-CAB-02**  
**11516-000023**

## Training Tools



**HeartSine samaritan PAD  
Trainer SAM 350P**

**TRN-350-XX\***  
**11516-000059**



**HeartSine samaritan PAD  
Trainer SAM 360P**

**TRN-360-XX\***  
**11516-000091**



**HeartSine samaritan PAD  
Trainer SAM 500P**

**TRN-500-XX\***  
**11516-000092**



**Trainer-Pak**  
(Replacement Electrode  
Cartridge)

**TRN-PAK-04**  
**11516-000017**

\* XX represents two-character country code.

# HeartSine Gateway™



**HeartSine Gateway**  
Includes removal tool, batteries  
and HeartSine Connected AED  
carrying case.  
**ACC-GTW-XX\***



**HeartSine Gateway  
Removal Tool**  
**ACC-REM-GW**



**HeartSine Gateway  
Batteries**  
**ACC-BAT-GW**

## Storage Options



**Wall Cabinet with Alarm,  
Custom Made**  
Dimensions:  
30.48 H x 27.94 W x 12.7 D cm  
(12 H x 11 W x 5 D in)  
Weight: 2.27 kg (5 lb)  
**PAD-CAB-04**  
**11516-000024**



**HeartSine AED Carrying  
Case**  
**PAD-BAG-01**  
**11516-000022**



**Mobile AED  
Rescue Backpack**  
**PAD-BAG-02**  
**11516-000114**



**Replacement electrodes**  
(for use with Trainer-Pak)  
**TRN-ACC-02** (10/pak)  
**11516-000009**  
**TRN-ACC-03** (25/pak)  
**11516-000011**



**HeartSine samaritan PAD  
Trainer Remote**  
**TRN-ACC-16** (SAM 350P)  
**11516-000014**  
**TRN-ACC-18** (SAM 360P)  
**11516-000016**



**TRN-ACC-17** (SAM 500P)  
**11516-000001**



**Battery Charger**  
(Replacement)  
**TRN-ACC-14**  
**11516-000012**

For further information contact us at [heartsinesupport@stryker.com](mailto:heartsinesupport@stryker.com) or visit our website at [heartsine.com](http://heartsine.com).



HeartSine Technologies, Ltd.  
203 Airport Road West  
Belfast, BT3 9ED  
United Kingdom  
Tel: +44 28 9093 9400  
Fax: +44 28 9093 9401



UL Classified. See complete marking on product.

The SAM 500P is not available for sale in the U.S.

All names herein are trademarks or registered trademarks of their respective owners.

© 2019 HeartSine Technologies. All rights reserved. H009-016-042-4 IE